Lansinoh.

HOT & COLD

POSTPARTUM THERAPY PACKS

Use cold for relief from pain and discomfort caused by labor and delivery.

Use hot for soothing relief from postpartum uterine cramping.

EN

TheraPearl® ColorSmart™ Technology turns the pearls white when heated or purple when chilled.* Follow heating or cooling instructions as directed. Apply pack for 20 minutes or as recommended by your physician.

For Cold Therapy:

Freeze pack for 4 hours before use.

For Heat Therapy - Microwave Instructions:

Microwave for the time indicated below. Gently massage pad to evenly distribute heat as indicated by color-changing pearls and test temperature prior to use. Add time in 5 second increments until desired temperature is reached. Heat only one pack at a time.

Suggested heating time for (1) one Postpartum Therapy Pack

25 seconds
at 700 watts

19 seconds
at 1000 watts

16 seconds
at 1250 watts

For Heat Therapy - Water Bath Instructions:

Fill a large pot with enough water to fully submerge the pack without touching the bottom of the pot. Bring water to a boil and remove pot from heating surface. **Do not immerse pack while water is being heated.** Use a blunt, nonmetallic tool to submerge pack into the water for 1-2 minutes. **Do not leave pack unattended while heating.**

Use a blunt, non-metallic tool to remove from water. If pack has not reached desired temperature, submerge again (in same water) for an additional minute.

Do not allow the pack to come in contact with the pot or heat source for prolonged time, excessive contact may rupture the pack. Towel dry pack before applying. Always inspect the pack for signs of damage or leaking prior to use and discard if damage is found.

*Follow heating and cooling instructions as directed. Always test the pack temperature before application and allow to cool if needed. Do not rely on color alone to indicate temperature.

Sleeve Instructions:

Place one pack either hot or cold into disposable sleeve. Disposable sleeves are intended for one time use. Packs may be used with or without the disposable sleeve.

Cleaning Instructions:

- After each use, thoroughly clean each pack with mild soap and water or 70% isopropyl alcohol. Rinse thoroughly until all soap and residue is removed.
- Allow to air dry completely prior to use.
- We do not recommend the use of ethanol, phenol, or chlorine-based products

Store in freezer or airtight container, out of direct light.

Suggested Uses

Perineum pain and hemorrhoids

Cold therapy, applied to the affected area can help reduce swelling and ease discomfort. Apply the pack for up to 20 minutes at a time allowing at least an hour in between uses or as recommended by your physician.

Postpartum uterine cramping

Heat therapy, applied to the affected area may help relieve the discomfort of cramping as your uterus begins to return to its original size. Apply the pack for up to 20 minutes at a time allowing at least an hour in between uses or as recommended by your physician.

CAUTION: Follow heating and cooling instructions as directed. Always test pack temperature prior to application.

Do not rely on color of beads to indicate temperature. Prolonged application may cause injury, including burns. Do not use on infants. Do not use if you have poor circulation. Do not use heat therapy on stitches or an open wound. Do not fall asleep using pack. Keep out of reach of children. For extreme, prolonged, or chronic pain, contact a physician. Use only as directed. Heating beyond suggested time may damage or destroy pack. Do not puncture and/or ingest contents. For external use only. Do not place disposable sleeves in microwave.

US Patents: <https://www.reckitt.com/patents>

Distributed by: / Distribuido por: Lansinoh Laboratories, Inc. Alexandria, VA 22314, USA
+1 800-292-4794

Lansinoh® is a registered trademark of Lansinoh Laboratories, Inc.
THERAPEARL® is a registered trademark of the Reckitt Group of Companies.
© 2022 Lansinoh Laboratories, Inc. All Rights Reserved. Made in China

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