Lansinoh®

SMARTPUMP 2.0

Double Electric Breast Pump

INSTRUCTIONS FOR USE

Intended Use of Breast Pump

The Lansinoh® Smartpump 2.0 Double Electric Breast Pump is intended to express and collect the breastmilk of a nursing woman for the purpose of feeding collected breastmilk to a baby. The Smartpump 2.0 is intended for a single user.

Indications of Breast Pump

Lansinoh Smartpump 2.0 Double Electric Breast Pumps are indicated to elimination of clinical situations such as milk stasis in the breast, engorgement and mastitis caused by insufficient milk discharge, relieve sore and cracked nipples and to bring out flat or inverted nipples. Also, Lansinoh Smartpump 2.0 Double Electric Breast Pumps enable mothers to provide breast milk for babies who cannot feed directly at the breast for a variety of reasons, such as latch-on problems, premature or unwell babies.

Contraindication of Breast Pump

There are no known contraindications.

Important Safeguards

Read All Instructions Before Using

Before first use, sanitize all parts that will come into contact with milk according to instruction on pg 10, section 3.1. Before each use visually inspect the individual components for cracks, chips, tears, discoloration, or deterioration. In the event that damage to the device is observed, please discontinue use until the parts have been replaced.

Danger To reduce the risk of electrocution:

- Always unplug this product immediately after using.
- Do not use while bathing.
- Do not place or store this product where it can fall or be pulled into a tub or sink.
- Do not place or drop pump into water or other liquid.
- Do not reach for an electrical product that has fallen into water.
Unplug immediately.

Warning:

The Lansinoh® Smartpump 2.0 Double Electric Breast Pump is a personal care item. To prevent cross contamination/infection, the Smartpump 2.0 tubing and breastmilk expression set should never be shared between mothers or re-sold from one mother to another. Each user should use her own Smartpump 2.0.

If you are a mother who is infected with Hepatitis B, Hepatitis C or Human Immunodeficiency Virus (HIV), pumping breastmilk using the Lansinoh® Smartpump 2.0 Double Electric Breast Pump will not reduce or remove the risk of transmitting the virus to your baby through your breastmilk.

When using electrical products, especially when children are present, basic safety precautions should always be maintained.

To reduce the risk of burns, electrocution, fire, or injury to persons:

- This product should never be left unattended when plugged in.
- Close supervision is necessary when this product is used near children.
- Use this product only for its intended use as described in this manual. Do not use attachments not recommended by the manufacturer.
- Repairs must be performed only by an authorized service agency. Do not repair yourself! No modifications to the device are permitted.
- Never operate this product if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water.
- Keep the cord away from heated surfaces.
- Never use while sleeping or drowsy.
- Never drop or insert any object into any opening or tubing.
- Do not use this product if it is broken or any parts are missing.
- Do not allow small children or pets to play with the pump unit, AC adapter, etc.
- Do not use outdoors with the AC adapter or operate where aerosol spray products are being used or where oxygen is being administered.

Caution:

If you feel pain in your breast or nipple during operation, turn the pump off or slide a finger between the flange and your breast to immediately break the suction. If irritation or discomfort occurs, discontinue use and see a doctor.

Save These Instructions along with your receipt for future reference.

Congratulations!

Your decision to provide your baby with breastmilk, even if you need to be apart, is one of the most precious gifts you can give your baby. The Lansinoh® Smartpump 2.0 is a safe, efficient, and effective choice for mom's who need to pump regularly.

Contents	Page
1. Quick Start Guide	5-6
2. Pump Diagram	7-9
2.1 Basic Controls, Display and Features	
2.2 Lansinoh® Smartpump 2.0 Breast Pump Parts	
3. Cleaning Instructions	10-11
3.1 Cleaning of Breastmilk Expression Set Before First Use	
3.2 Cleaning of Breastmilk Expression Set After Each Use	
3.3 Before Each Use of Breastmilk Expression Set	
4. Assembly	12-14
4.1 Assembling the Breastmilk Expression Set	
4.2 Connecting the Breastmilk Expression Set	
4.3 Connecting for Single Pumping	
5. Using the Lansinoh® Smartpump 2.0	15-18
5.1 Power Options	
5.2 Getting Started	
5.3 2-Phase Technology with Automatic Memory Function	
5.4 Customization for Comfort and Efficiency	
5.5 Completing your Pumping Session	
6. Pumping Tips	19
7. Troubleshooting	20
8. Breastmilk Storage Guidelines	21
9. Lansinoh® Bottle With NaturalWave® Nipple	22
10. Technical Data	23
11. Warranty	24-27

PLEASE DO NOT RETURN THIS PUMP TO THE RETAILER.

Lansinoh's dedicated Customer Service team is here to help.

If you experience any issues with your Lansinoh® Smartpump 2.0 breast pump please call **800-292-4794** or email **customerservice@lansinoh.com**, Monday - Friday 8am-8pm EST.

1. QUICK START GUIDE

Do NOT wash or sanitize tubing

Before first use, sanitize all parts that will come into contact with milk according to instruction on pg 10, section 3.1.

Please carefully read all instructions prior to use

Assembly

See detailed instruction on pgs 12-14. Proper assembly is essential to ensure proper function and suction.

1
Insert purple diaphragm with arrows pointing forward

2
Attach clear cover and twist so tubing connector faces back

3
Attach flange

4
Insert valve

5
Attach bottle

6
Attach tubing

7
Secure Y connector

8
Secure pump connector


9
Attach to pump base

Connect pump to power using AC Adapter or six AA Alkaline Batteries (not included)

Getting Started

Hold breastmilk expression set to your breasts so that the nipple is aligned to extend straight toward the nipple tunnel of the breast flange, ensuring there are no gaps for air to enter. Make sure set is tilted slightly downward to allow milk to flow into bottles.

2-Phase Technology Using Let-Down and Expression Phases

Press the power button  to start the pump in let-down phase. 

Let-down phase will run automatically for two minutes.

To begin expression phase sooner, press the let-down/expression toggle button  to move to the expression phase. Press the let-down/expression toggle button again at any point to switch back to let-down phase.

Customization for Comfort and Efficiency

Pumping Style Selection

Every baby feeds differently. The Smartpump 2.0 offers three different pumping styles to choose from, allowing you to choose the option that is most like your baby's feeding pattern. Change the pumping style using the style mode toggle button  to find the pumping style that works best for you.

Style 1

Style 2

Style 3

NOTE: When you turn the pump on for the very first time it starts in let-down mode, suction level 5. When you switch to the expression mode the pump will start in style 2, suction level 3.

Adjusting Suction Levels


Increase or decrease the suction to the maximum level that you find comfortable. Increase suction by pressing the  button. Decrease suction by pressing the  button.

NOTE: When increasing suction, the target suction level will flash, and the pump will slowly increase suction. The flashing will stop once the target level has been reached.


Lansinoh® is here to support you. If you have questions, please do not hesitate to **contact us at (800) 292-4794.**


2. PUMP DIAGRAM


2.1 Basic Controls, Display and Features

1. Power, On/Off button 
2. Suction Level Increase Control Button 
3. Pumping Style Toggle Button 
4. LCD Display
5. Let-Down/Expression Toggle Button 
6. Suction Level Decrease Control Button 
7. Let-Down/Expression Phase Indicator 
8. Time Display 
9. Pumping Style Indicator 
10. AC Power Indicator 
11. Battery Power Indicator 
12. Suction Level Indicator 


- 13. Tubing input
- 14. Battery compartment for 6 AA batteries (not included)
- 15. AC adapter input


AA Battery Compartment
(underneath pump)

AC Adapter

2.2 Lansinoh® Smartpump 2.0 Breast Pump Parts

Item	Item Number
A. Smartpump 2.0 Tubing	#53403
A1 Silicone Tubes 300 mm (2)	
A2 Silicone Tube 600 mm	
A3 Smartpump 2.0 Tubing Y Connector	#53405
A4 Smartpump 2.0 Pump Connector	#53453
B. Breastmilk Expression Set Includes:	
B1 ComfortFit® Breast Flange 25mm (2)	#50408
B2 ComfortFit® Breast Flange 30.5mm (2)	#53408
B3 Smartpump 2.0 Breast Flange Body (2)	#53454
B4 Smartpump 2.0 Diaphragm (2)	#53402
B5 Smartpump 2.0 Diaphragm Cap (2)	#53410
B6 White Valve (2 + 2 extra)	#53401
B7 5 oz (150 ml) Bottle (2)	#53451
B8 NaturalWave® Slow Flow Nipples (1)	#71027
C. Storage Lids (2)	#53450
D. Smartpump 2.0 AC Adapter/power cord	#50409
E. Wide Neck Bottle Holder (2)	#53452

Two extra white valves are included.
 To order replacement parts please visit www.Lansinoh.com

WARNING: To prevent cross contamination/infection, the Smartpump 2.0 tubing and breastmilk expression set should never be shared between mothers or re-sold from one mother to another. Each user should use her own Smartpump 2.0 tubing and breastmilk expression set.

3. CLEANING INSTRUCTIONS

Unless otherwise instructed by your healthcare provider, follow these sanitizing and cleaning instructions.

Do not boil or wash tubing or connectors!

WARNING: To prevent cross contamination/infection, the Smartpump 2.0 tubing and breastmilk expression set should never be shared between mothers or re-sold from one mother to another. Each user should use her own Smartpump 2.0.

When sanitizing parts, use enough water to keep parts afloat!

3.1 Cleaning of Breastmilk Expression Set Before First Use

1. Always wash and dry hands thoroughly with soap and water before handling the breastmilk expression set.
2. Disassemble the breastmilk expression set and place parts in a large pot.
3. Cover the parts in water making sure there is room for all of the parts to float freely in the pot.
4. Bring the water to a boil and allow the parts to boil for 5 minutes.
5. After boiling remove the parts from the pot with tongs and place on a clean towel to dry, covering the parts with a second towel.
6. Allow the parts to dry completely before assembly and use. Partially dried parts may affect the operation of the pump.
7. Once parts are dry they can be stored in a clean plastic bin with a lid or in a clean plastic zipper bag.

3.2 Cleaning of Breastmilk Expression Set After Each Use

1. Disassemble the breastmilk expression set, setting aside the tubing and connectors.
2. Rinse all parts with cold water immediately after each use. Place pump parts in a clean wash basin used only for washing infant feeding equipment. Do not place pump parts directly in the sink, because germs in sinks or drains could contaminate the pump.
3. Hand wash all parts of the breastmilk expression set using warm soapy water. A nipple brush may be used to scrub the hard plastic pieces of the breastmilk expression set, but should not be used on the soft silicone valve or diaphragm.
4. Rinse parts with running clean, warm water, or by submerging in fresh water in a separate basin that is used only for cleaning infant feeding items.
5. Place pump parts on a clean, unused dish towel or paper towel in an area protected from dirt and dust. Do not use a dish towel to rub or pat items dry. Allow the parts to dry completely before assembly and use. Partially dried parts may affect the operation of the pump.
6. Once parts are dry they can be stored in a clean plastic bin with a lid or in a clean plastic zipper bag.

7. Clean wash basin and bottle brush. Rinse them well and allow them to air-dry after each use. Wash them by hand or in a dishwasher at least every few days.

For Daily Sanitization of Breastmilk Expression Set:

1. Place clean parts in a pot and cover with water making sure there is room for all of the parts to float freely in the pot.
2. Bring the water to a boil and allow the parts to boil for 5 minutes.
3. After boiling remove the parts from the pot with tongs and place on a clean towel to dry, covering the parts with a second clean towel. Place wash basin and bottle brush on a clean, unused dish towel or paper towel in an area protected from dirt and dust. Do not use a dish towel to rub or pat items dry because doing so may transfer germs to the items.
4. Allow the parts to dry completely before assembly and use. Partially dried parts may affect the operation of the pump.
5. To use steam or microwave sterilizers: Follow the manufacturers instructions. Ensure all openings are face down in the sterilizer. Ensure you wash your hands before removing equipment from the sterilizer.

3.3 Before Each Use of Breastmilk Expression Set


1. Always wash and dry hands thoroughly with soap and water before handling the breastmilk expression set.
2. Before assembling the breastmilk expression set, be sure that all of the parts are completely dry.
3. Visually inspect the individual components for cracks, chips, tears, discoloration or deterioration. These may affect the functionality of the pump.
4. Inspect the tubing. If you find condensation in the air tube, wash and dry it immediately or replace the tubing. If you see breastmilk inside the tubing or tubing connectors do not attempt to wash, clean or sanitize the tubing or tubing connectors. Instead contact the manufacturer.
5. Should you have any questions or concerns about your parts please contact Lansinoh® at **800-292-4794** for further guidance.
6. To purchase replacement parts for your pump please visit **www.Lansinoh.com**

WARNING: It is highly recommended not to use any method other than the above mentioned cleaning methods.

WARNING: Use caution when cleaning the small white valve. It is important for proper suction. Do not use excessive force and do not clean with a baby bottle brush or a nipple brush. Gently wash by hand in hot, soapy water and rinse under hot running water.

4. ASSEMBLY


4.1 Assembling the Breastmilk Expression Set

1. Place purple silicone diaphragm, arrow side up, on top of breast flange body. Ensure diaphragm fits securely between rims of breast flange body and arrows point toward breast flange.


2. Place clear diaphragm cap over breast flange body and twist to secure so that the tubing portal is facing the back of the pump.

- 3.** Firmly push ComfortFit® Flange into breast flange body.

- 4.** Holding the base of the white valve, gently secure the valve into the bottom of the breast flange body.

NOTE: Use care when handling and cleaning the small white valve.

- 5.** Screw collection bottle onto breast flange body.

- 6.** Attach end of 300 mm tubing to back of diaphragm cap. Ensure that tubing is pushed all the way to the base of the attachment port and completely covers it.

- 7.** Attach other end of 300 mm tubing to tubing Y connector.

- 8.** Repeat assembly process for second breast flange body.

4.2 Connecting the Breastmilk Expression Set

Once Breastmilk Expression Set is assembled, connect to pump base:

1. Secure one end of 600 mm tubing to tubing Y connector.

2. Secure other end of 600 mm tubing firmly to the pump connector.

3. Insert pump connector into opening on left side of pump.

4.3 Connecting for Single Pumping

1. Remove one of the 300 mm tubing from the tubing Y connector.

2. Secure plug into open outlet of tubing Y connector.


5. USING THE LANSINOH® SMARTPUMP 2.0 DOUBLE ELECTRIC BREAST PUMP

5.1 Power Options

The pump may be operated using the following power supply modes:


AC Adapter

1. Connect the AC adapter to the power receptacle located on the right side of the pump. Plug the AC adapter into an available electrical outlet.
2. Electrical safety can only be achieved when the breast pump is connected to the AC power adapter supplied with the product. (See “Technical Data” section, page 23)

AA Alkaline Batteries (not included)

1. Make sure the batteries are inserted correctly according to the polarity (+) and (-) as indicated.
2. Close the compartment by gently pressing the cover until it snaps into place.


- Remove batteries if pump is not in use for a long period of time.
- Dispose of batteries according to regulations for environmental protection and recycling.
- Replace batteries as a simultaneous set. Do not mix old and new batteries.
- Do not mix alkaline and standard (carbon-zinc) batteries.
- If batteries are in the pump, but AC adapter is in use, the pump will use the adapter to power the pump, not the batteries.

For the most economical use, we recommend using your pump with the AC adapter provided. Your AC adapter is a universal AC plug adapter type and can be used with 100 Vac to 240 Vac systems with 50Hz and 60Hz mains.

5.2 Getting Started

If possible, wait until your breastmilk supply and breastfeeding are well established. This is normally at least 3-4 weeks after birth, unless otherwise advised by your healthcare provider.

1. Hold the breastmilk expression set to your breasts so that the nipple is aligned to extend straight toward the nipple tunnel of the breast flange.
2. Make sure the breast completely fills the breast flange so that a vacuum is formed and no air escapes.
3. The breastmilk expression set and nipple tunnel should tilt slightly downward to allow milk to flow naturally toward the bottle.

If you experience friction or discomfort on your nipples, you may need a different size flange. Larger size flanges are available. For more details, **please call 800-292-4794** or visit **www.Lansinoh.com**

5.3 2-Phase Technology with Automatic Memory Function

Using Let-Down and Expression Phases

1. Press and hold the on/off button to turn the pump on. The Lansinoh® Smartpump 2.0 will always start in the let-down phase.
NOTE: When you turn the pump on for the very first time it starts in let-down phase, suction level 5.
2. After two minutes, the let-down phase will automatically change to the expression phase.
NOTE: When you switch to the expression mode for the very first time the pump will start in pumping style 2, suction level 3.
3. If you prefer to bypass the pre-set two minute let-down phase or your milk lets down sooner than two minutes, press the let-down/expression toggle button to move to the expression phase.
4. Automatic memory function will save the last pumping style and suction level for each phase. However, the suction level will gradually build to the last used setting for a more comfortable pumping experience.

The style mode indicator will appear in the display to signify that you are in the expression phase.

5.4 Customization for Comfort and Efficiency

Using Customizable Pumping Styles

Choose Your Pumping Style

Each baby is unique – right down to the way he or she feeds. The Lansinoh® Smartpump 2.0 offers three different pumping styles to choose from – so you can select the one that most closely matches your baby’s feeding pattern.

	Pumping Style 1	
	Pumping Style 2	
	Pumping Style 3	

1. With the pump in expression mode, select the pumping style that is most comfortable and efficient for you by pressing the pumping style toggle button .
2. Pressing the button again will select the next pumping style.
NOTE: The pumping style toggle button only functions when the pump is in expression mode.
3. The mode will be indicated on the LCD display by the pumping style mode indicator.

As your baby grows, you may find that their feeding pattern changes. The customization features of the Smartpump 2.0 allow you to make adjustments so you're always getting the best results.

5.4 Customization for Comfort and Efficiency (cont.)


Tailor for comfort

Within each pumping style, there are eight levels of suction to choose from:

1. Adjust for stronger suction by pressing the suction level increase control button. (+)

Whenever the pump is switching to higher suction levels (within a mode or when switching modes) the targeted level will flash. The pump will slowly transition to the higher suction level. The flashing will stop once the target level is reached.

2. Adjust for lower suction by pressing the suction level decrease control button. (-)
3. The current suction level will be indicated by the number of dark bars on suction level indicator.

NOTE: If at any time the suction on your breast reaches uncomfortable levels, reduce the suction level by pressing the (-) button on the left side of the LCD screen or disconnect the tubing from its receptacle on the electric pump.

5.5 Completing Your Pumping Session

1. When a pumping session is finished, turn the pump off with the on/off button. (⏻)
2. After turning the pump off, the display will remain on for a few seconds.
3. Make sure there is no suction on the nipple(s) before removing the breast flanges from the breast. If necessary, break the suction by inserting a finger between the breast flange and the breast. Once the suction is broken, the flange can be removed safely.

PLEASE NOTE: The LCD display light will automatically switch off if no control buttons are pressed for 20 seconds in order to save power. If the pump unit is left on and unattended for more than 60 minutes, the control system will automatically power off. To turn unit back on, hold the power-on button for several seconds.

6. PUMPING TIPS

Pumping Tips

The most important part of successful pumping is the ability to trigger the let-down reflex. This is when your milk begins to flow freely. Warm, wet compresses applied to the breasts and gentle, circular massaging of the breasts prior to and/or during a pumping session can help encourage this natural reflex. Thinking of or looking at your baby, or a photo of your baby if you are apart, can also be helpful. Some moms can achieve more than one let-down during a pumping session. You may be able to achieve a second by returning to the let-down phase by pressing the Let-down/Expression Toggle Button.

How Long To Pump

Double pumping, compared to single pumping, will reduce the time you spend expressing your milk. Breastmilk is produced on a supply and demand basis. Expressing your breastmilk regularly will stimulate your body and help you maintain your milk supply. Just like with breastfeeding, it is not necessary to time your pumping sessions. As you get more experienced with pumping, you will get a better sense of the time it takes for your body to express your milk using a breast pump. If you have any breastfeeding or breast pump related questions, consult with an International Board Certified Lactation Consultant (IBCLC), breastfeeding counselor, or other healthcare provider such as your primary physician.

Pumping and Feeding Tips

- You may pump one breast while nursing the baby on the other. This allows you to maximize pumping effectiveness using the let-down reflex stimulated by your baby's sucking.
- Always test the temperature of any liquid before feeding your baby.
- Do not let a child take a bottle to bed or self-feed for long periods.
- To prevent possible choking hazard, test strength of the feeding bottle nipple by pulling on the bulb portion. Discontinue use if any tear or crack appears.

7. TROUBLESHOOTING

Troubleshooting the Lansinoh® Smartpump 2.0 Double Electric Breast Pump

Interruption of Electrical Service

The Smartpump 2.0 will shut off during a pumping session if the power source is changed while pumping.

This includes unplugging the power adapter or inserting the power adapter while pumping using batteries.

This safety feature is designed to prevent sudden changes in suction levels. If you lose power during a pumping session because you switched power sources mid-pumping session, simply press down on the power button to turn the pump back on.

If you experience an interruption in the normal operation of the electrical and/or mechanical function of the pump, simply do the following:

- Verify the power cord is still plugged in to the pump or the wall. If using batteries, verify all six batteries are fully inserted and none have come loose.
- If either the power cord or batteries became disconnected, fix the connection, and then press and hold the power button for two seconds to turn the pump back on.
- If the pump is shutting off and the adapter or batteries were firmly in place when this occurred, please contact Lansinoh Customer Service at **800-292-4794** or email **customerservice@lansinoh.com**.

Discomfort While Pumping

If you feel any discomfort, you may need a different size flange. Visit **www.Lansinoh.com** or call **800-292-4794**.

Decreased Suction

If it feels like the suction has decreased, try the following:

- Remove the tubing from pump and place a finger over the tubing connector.
- If you feel suction, the pump is operating properly, but the breastmilk expression set may not be correctly assembled or the white valve may be damaged.
- Make sure the breastmilk expression set has been properly assembled. Check the assembly information in the section “Assembling the Breastmilk Expression Set,” page 12-14.
- Inspect the white valve for a tear or pinhole. This component is critical to achieving proper suction. Replace the white valve with one of the spare valves provided with the product.

- If no suction is felt on your finger, try unplugging the adapter from the electrical outlet, wait a few seconds, and plug it back into the electrical outlet.
- When single pumping, check the tubing connector on the front of the pump to make sure it is closed for single pumping. If open, close the connector.
- If double pumping, ensure that the breastmilk expression set is securely connected on both breasts before starting the pump.

8. BREASTMILK STORAGE GUIDELINES

Breastmilk Storage and Thawing Guidelines for Full-term Healthy Babies

Check with your healthcare provider or an International Board Certified Lactation Consultant for specific storage instructions. When freezing, store in the back of the freezer and do not store milk in the door as temperature can fluctuate as the door is opened and closed.

Guide to storing fresh breastmilk for use with healthy, full-term babies

Place	Temperature	How Long	Things to Know
Countertop, table	Room temperature (up to 77°F)	Up to 4 hours	Containers should be covered and kept as cool as possible. Covering the container with a clean cool towel may keep milk cooler. Throw out any leftover milk within 2 hours after the baby is finished feeding.
Refrigerator	40°F	Up to 4 days	Store milk in the back of the refrigerator. When at work, it's OK to put breastmilk in a shared refrigerator. Be sure to label the container clearly.
Freezer	0°F or colder	Within 6 months is best. Up to 12 months is acceptable	Store milk toward the back of the freezer where the temperature is most constant. Milk stored at 0°F or colder is safe for longer durations, but the quality of the milk might not be as high.

Source: Adapted from 7th Edition American Academy of Pediatrics (AAP) Pediatric Nutrition Handbook (2014); 2nd Edition AAP/American College of Obstetricians and Gynecologists (ACOG) Breastfeeding Handbook for Physicians (2014); Academy of Breastfeeding Medicine (ABM) Clinical Protocol #8 Human Milk Storage Information for Home Use for Full-Term Infants(link is external) (2017); CDC Human Milk Storage Guidelines (2018).

Date all breastmilk and use the oldest breastmilk first.

Thawing

Thaw frozen milk in the refrigerator or place the containers under cool running water, gradually warming the water until the milk is thawed. Thawed milk can be kept in a refrigerator for up to 24 hours.

NOTE: DO NOT refreeze milk once it is thawed.

Breastmilk is not homogenized and you may notice that it separates as it sits. Gently swirl thawed milk before feeding to baby to mix the layers that have separated. Breastmilk can vary in color and consistency depending upon the time of the day that the milk was expressed, what you have eaten, and the age of the baby at the time of pumping.

NOTE: DO NOT MICROWAVE HUMAN MILK! Microwaving alters the composition of the milk. Microwaving also may create “hot spots” in the milk that could burn the infant.

NOTE: DO NOT thaw breastmilk in hot or boiling water.

9. LANSINOH® BOTTLE WITH NATURALWAVE® NIPPLE

Clinically Proven to Reduce Nipple Confusion in Breastfed Babies

The patented NaturalWave® Nipple enables baby to replicate natural sucking actions learned at the breast.

Exclusive breastfeeding is recommended for at least 3-4 weeks if possible before introducing a bottle and nipple.

Storage & Use

Prior to use, always inspect and pull the feeding nipple in all directions to check the product. Throw away at the first signs of damage or weakness. For safety and hygiene reasons, replace the nipple every 7 weeks.

Cleaning

Before first use, boil all parts for 5 minutes to sanitize. Ensure all parts are completely covered in water and free floating. Always wash all components of the product immediately after each use. The feeding nipple can be washed in hot soapy water and rinsed in clean water. Dishwasher safe on top rack only. Do not clean with solvents or harsh chemicals.

Assembly

Once thoroughly dry, insert feeding nipple through the top of the screw cap collar. Attach the screw cap collar securely onto the bottle neck. When carrying a full bottle, make sure it is closed with the collar in place, taking care to position the feeding nipple in the center of the collar to activate the “spill-proof” function of the product.

Precautions

Do not leave bottle in direct sunlight, heat or disinfectant (sterilizing solution), as this may weaken the nipple. Do not clean with or allow feeding nipple to come into contact with solvents or harsh chemicals. Damage could result. Fill with the required amount of liquid and tighten the screw cap before using. Do NOT overfill. Heating in a microwave is NOT recommended as it may produce uneven heating and high localized temperatures. If you choose to heat liquid this way, DO NOT SEAL BOTTLES. Swirl well before use, and always check the

Latching on

Baby's lips open outward and latch onto the areola



Peristaltic “wave-like” tongue movement

Baby's tongue moves in “wave-like” motion to compress nipple and extract milk

Swallowing

Back of baby's tongue rises, channeling the milk into the esophagus

temperature before feeding baby. Hot or unevenly heated liquid may burn baby's mouth. Do not warm liquid by placing the bottle directly onto gas or electric rings, warming plates, in electric ovens, etc. Wash after every use. Some liquids and sterile solutions may make the silicone become opaque. This does not alter any of its properties.

WARNING!

DO NOT MICROWAVE HUMAN MILK! Microwaving alters the composition of the milk and may also create "hotspots" in the milk that could burn the infant. Allowing your child to use this product for prolonged periods separate from regular mealtimes or to go to sleep while drinking liquids, other than water, from this bottle may cause or contribute to early childhood tooth decay. Always use this product with adult supervision. Keep all components not in use out of the reach of children. Do not leave baby alone with drinking equipment due to the risk of choking. Never use feeding nipple as a pacifier or soother. Always check liquid temperatures before feeding.

10. TECHNICAL DATA

1. Power Supply

Battery operation: 6 AA size alkaline batteries, all connected in series

Operating time: Minimum 30 minutes. Maximum 2 hours. Results are dependent on battery freshness and pump settings.

AC Adapter operation: 100-240VAC wall plug type AC adapter, galvanically isolated with 9Vdc output voltage

Pump power dissipation: 10.0 Watts maximum

Voltage regulation: Linear voltage regulator with low dropout voltage and short circuit protection

2. Control System

Control unit: Microcontroller with integral program memory

Program memory: 32K flash ROM, 2K RAM memory

HMI: LCD display with LED backlight, soft keypad

Single/double control: Single or double breast pumping possible

Selectable pumping style mode: 3 pumping style modes are selectable and controlled from keypad

Adjustable suction level: 8 suction levels are selectable and controlled from keypad

Suction vacuum level (Expression/Double pumping): Minimum 100 mmHg, Maximum 280mmHg.

Technical Data Bluetooth version 4.1

Range up to 10m

Power Requirements

AC Adapter (included)

AC Adapter – North American
Plug Configuration


100-240VAC wall plug type

AC adapter, galvanically isolated
adapter with 9Vdc output voltage.

Batteries (not included) 6 AA batteries,
1.5V alkaline

The expected use life of the pump is
500 hours.

**To maintain your power adapter follow
these instructions:**


DO NOT wrap cord of power adapter around the plug body.

CORRECT

DO NOT unplug power adapter by pulling on the cord.

CORRECT

DO NOT use power adapter if damaged.

11. WARRANTY

Limited Consumer Product Warranty (United States)

This product is warranted to the original consumer purchaser (the “purchaser”) to be free from defects in material and workmanship which are not commercially acceptable for the period of one year from the date of purchase on the pump mechanism (pump body and pump motor) and the period of 90 days from the date of purchase on all other product components. This warranty is provided only to products purchased in the United States and need parts shipped within the United States. Warranty coverage terminates if you (1) sell or otherwise transfer this product to another person; (2) improperly or incorrectly perform any maintenance or repair of this product; (3) use any non-Lansinoh® brand replacement parts with this product; or (4) misuse, abuse, or alter the product.

This warranty gives you specific legal rights and you may also have other rights, which vary, by location.

SYMBOLS	MEANING OF THE SYMBOLS
	Follow instructions for use
	Identifies the manufacturer
	Indicates the date of manufacture
	Type BF applied part
	Class II equipment
	Keep the device dry
	Keep the device away from sunlight
	Indicates the temperature limitation for operation, transport and storage.
	Indicates the humidity limitation for operation, transport and storage.
	Indicates the atmospheric pressure limitation for operation, transport and storage.
IP22	Indicates the protection against ingress of solid foreign objects and against harmful effects due to the ingress of water
	This pump is made of various metal and plastics. Before disposal, the device is to be rendered unusable and it must be disposed of as sorted waste in accordance with local waste electrical and electronic equipment regulations. To avoid the harmful effects on the environment and public health align with your defined local return and collection rules for waste electrical and electronic equipment (incl. batteries).
FC	Indicates the compliance with Federal Communications Commission requirements.
	Indicates that the packaging contains products intended to come in contact with food.
	Indicates the positioning of battery cell.
	Indicates to handle the fragile device with care.

LANSINOH® MAKES NO EXPRESS OR IMPLIED WARRANTY OF ANY KIND REGARDING THE PRODUCT OTHER THAN THOSE WARRANTIES EXPLICITLY SET FORTH HEREIN. ANY IMPLIED WARRANTY, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR ANY IMPLIED WARRANTY OF FITNESS FOR A PARTICULAR PURPOSE, TO THE EXTENT PERMITTED BY LAW, SHALL BE LIMITED IN DURATION TO A PERIOD OF 90 DAYS FROM THE DATE OF PURCHASE BY THE ORIGINAL PURCHASER.

In the event the product is found by Lansinoh® to not meet the above limited warranty, as purchaser's sole and exclusive remedy Lansinoh® will repair or, at the option of Lansinoh®, replace this product without charge for such replacement or parts or labor. No replacement or repair is provided outside of the United States.

NOTE: The purchaser shall bear all expense for returning this product to Lansinoh®.

Only Lansinoh® brand replacement parts are recommended for use with Lansinoh® breastfeeding products.

TO THE EXTENT PERMITTED BY LAW, LANSINOH® SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES INCLUDING, BUT NOT LIMITED TO, REPLACEMENT COSTS RESULTING FROM THE BREACH OF ANY EXPRESS OR IMPLIED WARRANTY.

If you wish to make a claim under this warranty, you must return the product to Lansinoh® Laboratories Inc., prepaid, together with your dated bill of sale or other dated proof of purchase and a brief statement of the problem, along with a Replacement Authorization (RA) number. You may obtain the appropriate Replacement Authorization number and shipping address by calling the phone number or sending an email to the email address provided below.

Items returned without an associated Replacement Authorization Number will be returned to sender.

To obtain the address to Return Goods or a Replacement Authorization Number:
800-292-4794 | customerservice@lansinoh.com

Declaration for Electromagnetic Compatibility (EMC)

For information use only

FCC Statement: This device complies with part 15 of the FCC rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency

energy, and if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Stacking: This device should not be used adjacent to or stacked with other equipment. And if adjacent or stacked use is necessary, this device should be observed to verify normal operation in the configuration in which it will be used.

Disposal: The symbol which can be found on the product or on the package surface indicates that the product shall not be disposed with other household wastes. It is the users responsibility to properly dipose of the product, or facilitate recycling activities. For further information about the recycling collection points, please contact your local authority or Lansinoh Customer Service.

EN 60601-1-2 STANDARD TEST REQUIREMENTS

Immunity Test Levels for ENCLOSURE PORT

Phenomenon	Basic EMC standard or test method	IMMUNITY TEST LEVELS	
		Professional healthcare facility environment	HOME HEALTHCARE ENVIRONMENT
ELECTROSTATIC DISCHARGE	IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	
Radiated RF EM fields ^{a)}	IEC 61000-4-3	3 V/m ^{d)} 80 MHz - 2,7 GHz ^{b)} 80 % AM at 1 kHz ^{c)}	10 V/m ^{d)} 80 MHz - 2,7 GHz ^{b)} 80 % AM at 1 kHz ^{c)}
Proximity fields from RF wireless communications equipment	IEC 61000-4-3	See 8.10	
RATED power frequency magnetic fields ^{a) e)}	IEC 61000-4-8	30 A/m ^{g)} 50 Hz or 60 Hz	

Immunity Test Level for Input A.C. Power PORT

Phenomenon	Basic EMC standard or test method	IMMUNITY TEST LEVELS	
		Professional healthcare facility environment	HOME HEALTHCARE ENVIRONMENT
Electrical fast transients/bursts ^{a) 1) c)}	IEC 61000-4-4	± 2 kV 100 kHz repetition frequency	
Surges ^{a) b) 1) c)} Line-to-line	IEC 61000-4-5	± 0,5 kV, ± 1 kV	
Surges ^{a) b) 1) k) c)} Line-to-ground	IEC 61000-4-5	± 0,5 kV, ± 1 kV, ± 2 kV	
Conducted disturbances induced by RF fields ^{c) d) e)}	IEC 61000-4-6	3 V ^{m)} 0, 15 MHz – 80 MHz 6 V ^{m)} in ISM bands between 0,15 MHz and 80 MHz ⁿ⁾ 80 % AM at 1 kHz ^{e)}	3 V ^{m)} 0, 15 MHz – 80 MHz 6 V ^{m)} in ISM and amateur radio bands between 0,15 MHz and 80 MHz ⁿ⁾ 80 % AM at 1 kHz ^{e)}
Voltage dips ^{f) a) 1)}	IEC 61000-4-11	0 % UT ; 0,5 cycle ^{q)} At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° ^{q)}	
		0 % UT ; 1 cycle and 70 % UT ; 25/30 cycles ^{h)} Single phase: at 0°	
Voltage Interruptions ^{f) 1) c) 1)}	IEC 61000-4-11	0 % UT ; 250/300 cycles ^{h)}	

Conducted Emission Test Level

Frequency range MHz	Quasi-peak dB(µV/m)	Average dB(µV/m)
0,15 – 0,50	66 Decreasing linearly with logarithm of frequency to 56	56 Decreasing linearly with logarithm of frequency to 46
0,50 – 5	56	46
5 – 30	60	50
At the transition frequency, the more stringent limit shall apply.		

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Manufacturer: Lansinoh® Laboratories Saglik Gerecleri Tasarim San. Tic. Ltd. Sti.